



CCP5 Summer School

July 2023

Practical sessions: Structure-based and ligand-based design of C5 inhibitors with BioSolveIT suite

DR. AGNIESZKA K. BRONOWSKA

SCHOOL OF NATURAL ENVIRONMENTAL SCIENCES

NEWCASTLE UNIVERSITY

Contents

- Molecular docking and virtual screening
- Introduction to SeeSAR
- Application 1: lead optimisation
- Application 2: core expansion
- Application 3: core replacement/scaffold-hopping
- Application 4: fast generation and evaluation of analogues
- Application 5: Virtual screening of massive virtual libraries with InfiniSee (Scaffold Hopper) and SeeSAR

Molecular docking: two main tasks

 Sampling of ligand conformational space and pose generation (geometry)

Scoring protein-ligand complexes (energetics)


Molecular docking: to-do list

Problem: a pair of molecules represented by their 3D coordinates

- Decide whether the molecules will form a complex;
- Determine the binding affinity (free binding energy);
- Predict the 3D structure of the complex (binding mode);
- Deduce function (agonist/antagonist)*;


Installing SeeSAR and InfiniSee

- Platform-independent
- Quick to install
- Temporary license key (works for both) provided





Download SeeSAR for your operating system (System requirements):



Introduction to SeeSAR

Tasks 1 and 2

•Import your target structure to SeeSAR


You can either import directly from RCSB PDB Data Bank, or read your local file


Assign the binding site

The process is automated: user defines the site either by existing ligand (docked or experimental), or by unoccupied pockets

For the covalent docking, you need to prepare your site (target residue) – or to start with the ligand covalently bound

Getting started






The binding site definition in SeeSAR (1)


PDB code: 1QCF

1.


You must tick here ————

2.






3. The outcome




The binding site definition in SeeSAR (2)

Using unoccupied pockets







Lead optimisation

PDB code (source): 8AYH

Use PDB file provided: C5_H1H_CCP5.pdb

CryoEM structure (3.35 Å resolution) of human complement C5 in complex with small molecule inhibitor (H1H) and CVF


 IC_{50} : 0.1 - 5 nM, from 3 assays

- •Task 3: open the target-inhibitor complex in SeeSAR, indicate the binding site by the inhibitor (H1H), and bring it to the Molecule Editor mode
- Once you are in the Editor mode, you can manually modify/grow your compound
- •Before you start modifying: what is the affinity range calculated? How does this compare to the experimental binding affinity?
- What you may try (a good practice): molecular docking
- Bring your compound to the **Docking mode**


Complement C5 inhibitor: docking

Task 4: Use standard (noncovalent) docking, 50 poses, high clash tolerance and maximum allowed ring conformations

First you generate the poses, then you select them (select all) and rank them by calculating their binding affinities



Task 5: Sort the poses from the best to the worst and bring the best one back to the Editor. Compare it with the starting structure. What do you see?



Task 6: Try 5-10 modifications to improve the binding affinity (5-10 minutes).

Do not forget to save the project file every now and then.

Task 7: Bring an "external" molecule to SeeSAR (Editor mode) and dock it to C5.


(I will let you to figure out how)

Remember: SeeSAR can deal with a lot of file formats – including SMILES strings.

Good practice: after scoring and ranking, delete the worst/clashed poses to avoid massive files


Core expansion in SeeSAR


Configuring Inspirator mode



Task 8: Core expansion in SeeSAR


- •We will select a central core of H1H for the expansion
- •Bring H1H compound from the Editor to the Inspirator mode








Inspirator







5	52	☆	₪	4	@	Name	Estimated Affinity						Fragment	LLE	Tor.	Int cla
	~	^	~	•	,		pM	nM		μМ	mM		Origin			Cla
1						H1H_A_768	\vdash	4				n.a.		8	•	
2						H1H_A1324_5						0.00	CSSB00011324			
3						H1H_A6515_5						0.00	CSSB00186515			
4						H1H_A7244_5						0.00	CSSB00057244			
5						H1H_A1346_5						0.00	CSSB00741346			
6						H1H_A0718_5						0.00	CSSB00560718			
7						H1H_A9724_5						0.00	CSSB00729724			
8						H1H_A9698_5						0.00	CSSB00989698			
9						H1H_A4871_5						0.00	CSSB00724871			
10						H1H_A2727_5						0.00	CSSB00582727			
11						H1H_A0883_5						0.00	CSSB00560883			

•5				Name Estimated Affinity				Deviation Fragment		LLE	Tor.	Int					
		ν.	ш	7	a)	3		pM	nM		μМ	mM		Origin			cla
1	\checkmark						H1H_A6515_5	一					0.00	CSSB00186515			
2							H1H_A_768	_					n.a.		8	•	
3							H1H_A1324_5	H					0.00	CSSB00011324	(B)		
4							H1H_A4871_5	Н	—				0.00	CSSB00724871	•		
5							H1H_A0883_5	F	_				0.00	CSSB00560883	(%)		
6							H1H_A9698_5	F	—				0.00	CSSB00989698	9 }		
7							H1H_A0718_5	- 1	$oldsymbol{}$				0.00	CSSB00560718	9 }		
8							H1H_A1346_5		_	ı			0.00	CSSB00741346	(4)		
9	$\overline{\mathbf{Z}}$						H1H_A9724_5		\vdash	•			0.00	CSSB00729724	•		
10							H1H_A7244_5		—	-			0.00	CSSB00057244			
11							H1H_A2727_5		—	–			0.00	CSSB00582727			

Inspirator: example results





Inspirator

It may be useful to output the

following:


•MW

•logP

TPSA

Torsional quality estimates

Clashes




Task 9: Core expansion in SeeSAR


- Repeat the expansion using different direction of the growth
- •You may start from H1H, of your "best" analogue from the previous growth iteration

Core replacement

Task 10: Core replacement in Inspirator

- Requires 2 bonds selected
- Does not cut through the ring
- May be done iteratively with core expansion/growth









Generation of analogues

SEESAR VERSION 13+ ONLY

Task 11: Analogues in Inspirator

•Requires a binding site and at least one bound molecule





Virtual Screening

Virtual screening using SeeSAR and InfiniSee

- •Docking of a large virtual libraries (1,000+) of compounds
- •Libraries: collections of small molecules for virtual screening
- •Sources: open-access (e.g. ZINC, ChEMBL) or commercial (e.g. Enamine) virtual libraries
- Types of libraries commonly used: diversity sets, target-focused, custom

Building your own C5 focused library in InfiniSee

- You may use H1H as a query ("bait")
- Load it into InfiniSee (Scaffold Hopper mode)
- Choose the chemical space




- You can pull thousands of synthetically-feasible analogues in minutes to create your own bespoke focused libraries (SDF/SMILES formats)
- Next, you can virtually screen these analogues


Scaffold Hopper in InfiniSee




Configuring InfiniSee




Customise similarity features and results








InfiniSee/Scaffold Hopper: results





Chemical space search and virtual screening

Task 12: Using you "best" C5 inhibitor as a "bait" (query), find 100 compounds in GalaXi space.

Out of these 100, pick 10, save them as either SDF or SMILES, upload to SeeSAR and virtually screen them against C5.

Of all compounds virtually tested for C5 binding: select your favourite, export the complex as the PDB file, and you may follow it up using all-atom MD simulations.

Recap: limitations and known issues

- •Protein is considered rigid: ideally, you should follow your calculations by running short MD simulations on the complexes and recalculating the binding affinities
- •Results are very sensitive to even small changes in the conformation of the protein
- •Binding affinities for certain groups are not reproduced well: hydrophobic effect tends to be overestimated, while highly polar groups are underestimated
- •Every now and then, weird protonation states suggested (you can always manually adjust) and med-chem nonsense molecules suggested in core expansion in Inspirator
- Workflows are limited to small molecules
- Med-chem properties and/or synthetic feasibility of suggested analogues may be problematic
- •Technical: 50,000 compounds/rows in SeeSAR GUI (you need to use KNIME to "downsize" very large data sets, or use non-GUI version)